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DESCRIPTION

It is a monofilament consisting of a sterile crystalline isomer.
The linear hydrocarbon Polymer does not contain unsaturated bonds.
The suture is colored with D&C#6 blue OR copper phthalocyanine blue.

INDICATION

The Monofilament Polypropylene surgical suture can be used in all tissues recommended by the USP pharmacopoeia for the use of non-absorbable sutures. Due to its relative biological inertia, it is recommended for use in tissues where minimal tissue reaction is desired. Because it is non-capillary and resistant to infection development, it can be used in infected or contaminated tissues to minimize or eliminate delayed sinus drainage and suture extrusion.

AVAILABILITY

Available in USP diameters from 11-0 to 2. Supplied sterile in pre-set sizes, with a swaged needle, blue-colored. Boxes of 6, 12, 24 and 36 units.

INSTRUCTIONS FOR USE

The Monofilament Polypropylene surgical suture should be used according to the procedure.
The size and diameter of the suture thread and needle to be used in the procedures will be established at the surgeon's discretion. Apply the suture on the edges of the surgical incision, using the "interrupted" or "continuous" techniques according to the surgical techniques of incision closure and medical criteria.
The Monofilament Polypropylene surgical suture is not absorbed and must, therefore, be removed (short term use).

PRECAUTIONS, WARNINGS, SPECIAL CARE AND CLARIFICATIONS ON THE USE OF THE PRODUCT

The product should only be used by or under the direction of a surgeon, according to applicable surgical techniques and should only be handled by qualified personnel. It is recommended that proper surgical practices be followed for the care of infected or contaminated wounds.

Care should be taken in the handling of this product to ensure that it is done aseptically. Care should be taken in the handling of surgical needles to avoid accidental puncture.

Like any other suture, care must be taken in its handling, avoiding damage, especially for being a monofilament. Kneading, scraping or binding of the suture with rough, sharp or cutting instruments must be avoided. Do not hold the suture with unprotected needle driver and needle holders. An exception is made when knotting with instruments, when the thread must be held at its extremities.

The guarantee of appropriate knots is given by the surgical circumstances and by the surgeon's experience, considering surgical techniques of flat knots followed by a square knot with several superposed ones.

Dispose of used needles in proper recipients.

Disposal of used sutures: danger of biological contamination. Adopt disposal measures, considering the risk of biological contamination and perforating material.

Disposal of unused sutures: adopt disposal measures for piercing-cutting material.

Do not use the product for purposes other than those recommended by the manufacturer.

Do not use the product after its expiration date.

Do not resterilize. Discard opened and unused product; discard leftover suture.

Do not use if package is violated or wet.

Reprocessing forbidden. Single-use health product.

CONTRAINDICATIONS

The Monofilament Polypropylene surgical suture is contraindicated for use in ophthalmic, neurological, cardiovascular and microsurgery surgeries.

ADVERSE REACTIONS

Adverse reactions reported from the use of this product include wound dehiscence, minimal inflammatory tissue reaction and temporary irritation at the wound site.

CONSERVATION

Transport and store at room temperature (15 to 30°C), in a cool, dry place, protected from humidity. Keep away from direct sunlight, heat sources, dust and bad weather conditions.

SINGLE USE PRODUCT. DESTROY IT AFTER USE.

STERILE PRODUCT - ETHYLENE OXIDE.

Validity 5 years from date of manufacture.

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